

LorX[®]

tlif peek cage



 **Tria Spine[®]**
SPINE | ORTHOPEDICS

DESIGN FEATURES

Better Stability

enhanced instrumentation for secure insertion



Special Teeth Design

specially designed teeth structure for easy insertion and secure placement




Titanium Pins

three titanium markers for the verification of anterior and posterior implant placement

Biocompatibility

radiolucent and biocompatible PEEK material (VESTAKEEP by Evonik, ASTM F2026) which provides modulus of elasticity similar to the bone



Enhanced Bone Fusion

two adequate graft spaces for bone fusion

Anatomical Shape

original and anatomical design

LorX[®] TLIF PEEK Cage

LX-27-317	11 x 31 mm	7 mm	LX-27-3110	11 x 31 mm	10 mm
LX-27-318	11 x 31 mm	8 mm	LX-27-3111	11 x 31 mm	11 mm
LX-27-319	11 x 31 mm	9 mm	LX-27-3112	11 x 31 mm	12 mm